


# Getting Started with EPICS

## *Applications / Special Topics*

### Introduction to synApps (v5.1)

*Tim Mooney*  
1/11/2005

## **Argonne National Laboratory**



A U.S. Department of Energy  
Office of Science Laboratory  
Operated by The University of Chicago



# What is synApps?

---

- A collection of EPICS applications for synchrotron-beamline users  
<http://www.aps.anl.gov/aod/bcda/synApps>
  
- **EPICS modules and build/configuration tools:**
  - Modules: *autosave, calc, camac, ccd, dac128V, dxp, ip, ip330, ipUnidig, love, mca, motor, optics, quadEM, sscan, std, vme, xxx*
  - Build/config: *config and utils directories*
  
- **Related clients, libraries, and visualization tools:**
  - IDL: *scanSee, mca display, ezcaIDL, ezcaScan, ez\_fit, HDF translator/browser, Ascii-format plotter, image processors, etc.*
  - CA-Server based CCD control
  - some python support

# synApps modules

---

- **Modules contain the following kinds of support:**
  - Compiled code; libraries
    - *E.g., record and device support*
    - *State-Notation-Language programs*
  - EPICS databases and autosave-request files
    - *A database is a **program** written in a high-level language.*
    - *One or more copies of a database can be run, each with its own private variables (PV's).*
    - *The database designer recommends PV's to be autosaved by naming them in a .req file; you can override with a private copy of the file.*
  - MEDM-display files
    - *The default user interface*
  - Documentation

# Other EPICS modules used by synApps

---

- **asyn 4.1**
- **ipac 2.8**
- **seq 2.0.8 (9?)**
- **genSub 1.6**
- **vxStats 1.7.2c**
- **allenBradley 2.1**

# autosave module

---

- **Records latest values of selected EPICS PVs; restores those values when the ioc restarts.**
  - not an archiver; only the latest value is saved
  - not the same as saveData, which writes scan data
  - When a list of PV's is saved, the entire list is written, even if only one PV has changed.
- **Can save/restore any scalar or array-valued PV (synApps 5.1)**
  - Array-valued PV must be hosted by the ioc that does the restore operation. (Typically, all ioc's save/restore their own PV's.)
  - DBF\_MENU, DBF\_ENUM PV's are handled by number.
- **Save operation uses channel access for scalars.**
- **Restore operation uses static database access for scalars.**
- **Arrays are saved and restored with database access.**

## ...autosave module

---

- **Three restore options for save files:**
  - 1) before record/device initialization
    - Motor positions must be restored at this time.
    - Arrays cannot be restored at this time. \*
    - PV's that are DBF\_NOACCESS before record init (e.g., genSub variable-type fields) cannot be restored at this time. \*
  - 2) after record/device initialization
    - to override record-initialization values
    - Link fields cannot be restored at this time. \*
  - 3) both before and after record initialization
    - The 'auto\_settings.sav' file is restored at both times.
    - It's not an error to attempt to restore a PV at the wrong time.
    - If you restore a motor position at this time, you override the value read from hardware, without writing to hardware.

\* Not illegal, just doesn't work

## ...autosave module

---

- **PV lists can use include files** (e.g., `<database_name>.req`), include path.
  - Database developer can supply default include file with database.
  - User can override with custom include file.
- **Save triggers:**
  - on change of any PV in the list
  - periodically
  - on change of a trigger PV
  - manual
- **User can reload save sets.**
- **Autosave can recover from file-server reboot (synApps 4.6+).**
  - Currently, only on vxWorks
- **User can choose to save redundant files (synApps 4.6+).**
- **Autosave reports status via EPICS PV's (synApps 5.1+).**

# ...autosave module

- Sample request file

```

xxx:my_PV.VAL
xxx:my_array_PV.VAL
file motor_settings.req P=$(P),M=m1
...
<END><lf>

```

*PV names*  
*Keyword*  
*Macro replacements*  
*Name of include file*

- Sample save file

```

# save/restore V4.4 Automatically generated...
xxx:my_PV.VAL 1.0
xxx:my_array_PV.VAL @array@ { "0" "0.1" ... "10.2" }
xxx:m1.DIR 0
xxx:m1.DHLM 100
xxx:m1.DLLM -100
...

```




# calc module

---

- **Support for evaluation of string or numeric expressions entered at run time (or at database-configure time)**
- **Records**
  - **sCalcout** – like calcout, but also supports string expressions; user can specify wait-for-completion.
  - **swait** – like calcout, but uses recDynLink (no “PP MS” link attributes)
  - **transform** – like 16 calcout records that share a PV data pool
- **Other code**
  - string-calc engine
  - sCalcout soft device support (with wait-for-completion option)
  - interpolation routines for genSub record
  - (yet another) averaging routine for sub record

# ...calc module

- Databases, medm displays for run-time programming
  - userCalc
  - userStringCalc
  - userTransform
  - userAve
  - arrayTest
  - interpolation
- Examples of **ALL** calc expressions (normal and stringCalc) can be found in synApps MEDM help displays



The image shows two overlapping windows from the EPICS calc module help displays. The background window is titled 'calcExamples.adl' and lists various mathematical and string functions. The foreground window is titled 'calcMiscExamples.adl' and shows examples of more complex expressions.

**calcExamples.adl**

**Algebraic functions/operators**

ABS(a)	LOG(a)	SQR(a)
CEIL(a)	LOGE(a)	SQRT(a)
EXP(a)	MAX(a,b,...)	a^b
FLOOR(a)	MIN(a,b,...)	a**b
INT(a)	NINT(a)	a+b
LN(a)	NOT(a)	a-b

**Trigonometric functions/operators**

ACOS(a)	COS(a)	SINH(a)
ATAN(a)	COSH(a)	TAN(a)
ATAN2(a,b)	SIN(a)	TANH(a)

**Relational operators**

a>=b	a<=b	a!=b
a>b	a<b	a#b

**Bitwise operators**

a b	a&b	a XOR b
a OR b	a AND b	~a

**Misc operators & symbols**

(a)	a?b:c	@a
RNDM		

**String functions/operators**

BYTE(aa)	PRINTF('%f',aa)	aa+bb
DBL(aa)	SSCANF(aa,'%lf')	aa>=bb
INT(aa)	\$(aa,'%lf')	aa>?bb
MAX(aa,bb,...)	STR(aa)	aa!=bb
MIN(aa,bb,...)	aa[1,3]	aa>bb
NINT(aa)	aa[bb,cc]	aa==bb
\$P('%f',aa)	'abcdef'	aa- bb

NOTE: can use single or double quotes at but only single quotes in a static database  
'a', 'b', etc. are numeric arguments.  
'aa', 'bb', etc. are string arguments.  
Items in blue are useable only in transform and sCalcout records.

**calcMiscExamples.adl**

TOP ALG TRIG REL BIT MISC STR

**Miscellaneous expressions**

**Grouping**

(a)  
CAN'T USE '[' or '{'

(A)ANDB  
SAME AS A&B. BUT NOTE...

AAANDB  
SCALCOUT RECORD WILL PARSE 'AA' AS A STRING VARIABLE

**'?:' (if-else) operator**

A?B:C  
IF (A!=0) RESULT IS B, ELSE C

**Array expressions**

@0  
GET VARIABLE 0 (<@0=A, @1=B, ...)

@A+3  
GET VARIABLE A; ADD 3 TO IT

@(A+3)  
GET VARIABLE A+3

@@0  
GET STRING VAR 0 (<@@0=AA, ...)

**Random number**

RNDM  
RANDOM VALUE IN RANGE [0,1]

# camac module

---

- **Communication with CAMAC crate/modules**
- **Records**
  - camac – generic BCNAF/data for run-time camac control
- **Devices supported**
  - VME bus adapter
  - CAMAC crate controller
  - E500 motor controller
  - RTC-018 real-time clock
  - QS-450 quad scaler
  - DXP spectroscopy system (now in *dxp* module)

# ccd module

---

- **Support for area detectors (CCD's and image plates)**
- **Supported devices**
  - MAR 165 CCD
  - MAR 345 image-plate reader
  - Roper (all WinView-supported CCD's, including former Princeton and most former Photometrics devices)
  - Bruker SMART CCD
- **Can control, at minimum**
  - exposure time
  - file name
  - data-acquisition start
  - wait for acquisition to complete
  - much more for most devices
- **See lecture “*Detectors and Feedback.*”**

# config directory

---

- **Configures and builds all modules in or used by synApps**
- **MASTER\_RELEASE**
  - specifies version number and file path to EPICS base, and to every module in or used by synApps
- **makeReleaseConsistent.pl**
  - Edits <module>/configure/RELEASE for every module in or used by synApps, to agree with MASTER\_RELEASE
  - “gnumake release” causes this to run.
- **Makefile**
  - “gnumake <whatever>”, in config directory, does <whatever> for all modules.

# dac128V module

---

- **device support, database, and MEDM displays for dac128V IndustryPack module**
  - 8-channel, 12-bit DAC
  - Support exists to run a DAC channel manually, or according to an algorithm written at run time, or as a *scan* positioner, or as part of a PID feedback loop.
  
- **See lecture “*Detectors and Feedback.*”**

# documentation directory

---

- ***TOP-level synApps documentation***
  - What synApps is
  - How to build it
  - How to make a user application from the 'xxx' sample module
  - How to fit the user application to a particular set of hardware
  
- **This presentation**

# dxp module

---

- record, device support, databases, and MEDM displays for XIA DXP and Saturn spectroscopy systems
- dxp record for setting DXP parameters
- device support for the mca record
- See lecture “*Detectors and Feedback.*”



# ip module

---

- **device support, SNL code, databases, and MEDM displays for *many* message-based devices**
  - originally, for devices supported via IndustryPack hardware
  - Note some of this support will inevitably be out of date -- pending access to hardware for testing.
  
- **deviceCmdReply (was serial\_OI\_block, GPIB\_OI\_block)**
  - Used to write support at run time for one command/reply message
  - sCalcout to format output string
  - asyn record to write/read device
  - sCalcout record to parse reply
  
- **devXxStrParm device support**
  - probably will be replaced by streams/asyn

# ip330 module

---

- device support, databases, and MEDM displays for the IP330 ADC IndustryPack module
- 16/32 channel, 16-bit ADC
  - ip330Scan for periodic, averaged reads of ADC channels
  - ip330Sweep, with the MCA record, for using ip330 as a waveform-digitizer
  - ip330PID for using the ip330 in a fast-feedback loop
- See lecture “*Detectors and Feedback.*”

# ipUnidig module

---

- device support, databases, and MEDM displays for the IPUnidig digital I/O IndustryPack module
- IP-UD-I 24-channel input/output/interrupt module
- DIO316I 48-bit digital I/O module
- See lecture “*Detectors and Feedback.*”

# love module

---

- **Support for Love controllers**
  - orphaned, currently under re-development for EPICS 3.14

# mca module

---

- **Support for multichannel analyzers, multichannel scalars, and other array-valued detectors**
- **mca record**
- **device support**
  - Canberra 556 AIM module (MCA and ICB controller)
  - DSA-2000 Ethernet MCA
  - various Canberra-ICB modules for spectroscopy
  - SIS 3801 (Struck STR7201) MCS
  - (DXP support in dxp module)
  - (IP330 support in ip330 module)
  - (quadEM support in quadEM module)
- **See lecture “*Detectors and Feedback.*”**

# motor module

---

- **Motor record and device support**
  - stepper and servo motors
  - soft-motor support
    - *Put motor “face” on, e.g., a DAC channel*
    - *Drive a hard motor through a nonlinear transform*
  - user/dial/raw coordinates
  - backlash-takeout algorithm
  - pre/post move commands
  - many more features
  
- **See lecture “*Motors.*”**

# optics module

---

- **Slits and mirrors**
  - *Four virtual positioners; two real motors*
  - *Automatic sync to motor positions*
  - *Completion reporting*
- **Monochromators**
  - *Nondispersive double-crystal*
    - Geometries: (Y1, Z2), (Y2, Z2)
    - Crystal species: Si, Ge, Diamond, Si (77K)
    - Miller indices, allowed reflections
    - Operational modes:
      - *Use/Set*
      - *Manual/Auto*
    - Managing the vertical beam offset
    - Automatic sync to motor positions

# ...optics module

---

- **...Monochromators**
  - *Spherical grating*
    - Geometrical variables:
      - 1) *Grating line density; radius*
      - 2) *Tangent-arm length*
      - 3) *Diffraction order*
      - 4) *Input/output slit distances*
    - Operational modes:
      - *Use/Set*
      - *Manual/Auto*
    - *Grating-stripe list*
    - *Manual sync to motor positions*



# ...optics module

---

- **...Monochromators**
  - *Dispersive double-crystal*
    - Geometries: nested, symmetric
    - Crystal species: Si, Ge, Diamond, Si (77K)
    - Miller indices, allowed reflections
    - Operational modes:
      - *Use/Set*
      - *Manual/Auto*
      - *Theta1 / Theta1&2 / Rock Theta2*
    - Accommodate incident-beam angle shift (“world offset”)
    - Automatic sync to motor positions

# ...optics module


---

- **Optical table**
  - **Table** record supports a six-degree-of-freedom optical table.
  - User/client can write either to  $(x, y, z, \theta_x, \theta_y, \theta_z)$ , or to underlying motor records.
  - Table rotates about user-specified point.
  - Table database includes a list of rotation points, selected by menu.
  - Can recover table position from motor positions
  - Partial support for fewer than six degrees of freedom

# ...optics module

- **Optical table**

- Geometries
  - SRI
  - GeoCARS
  - Newport
  - PNC



- Calibration/sync

- Use/Set – changes to [X, Y, ..]. move table / change calibration
  - Zero – redefine current [X, Y, ...] as zero
  - Sync – update [X, Y, ...] from motors, honoring calibration
  - Init – clear calibration and sync to motors
- 
- Table record sets motor speeds so that motors start/stop together.

# quadEM module

---

- Support for APS Detector Group's (Steve Ross) four-input electrometer.
- See lecture “*Detectors and Feedback.*”

# sscan module

---

- **Support for user-programmable data-acquisition**
  - **sscan** and **busy** records
  - saveData
  - recDynLink
- **A one-dimensional scan:**
  - Do NPTS times:
    - Set conditions e.g., move motors; wait for completion
    - Trigger detectors e.g., start scaler; wait for completion
    - Acquire data read detector signals; store in arrays
  - Write data to NFS file
- **Multidimensional scan:**
  - Same as a 1-D, but detector trigger executes inner-loop scan.
  - saveData monitors a set of **sscan** records, determines scan dimension when scan starts, and writes data as it is acquired.

# ...sscan module

---

- **scan features:**
  - Three 1-D scan types: constant-step-size, table-driven, fly
  - Unlimited number of data points, scan dimensions
  - 0-4 positioners, 0-4 detector triggers, 0-70 detector signals
  - Acquisition from scalar and 1-D-array-valued PV's
  - Detector/client wait, data-storage wait
  - Pause/resume, abort
  - Double buffered: can write 1-D acquired data during next 1-D scan
  - *saveData* writes self-describing XDR-format (".mda") files to NFS-mounted disk (vxWorks only, at present).
  - A positioner can have private scan parameters (scanparm record).
  - After-scan actions include move to peak, valley, and edge.
  - scanparm record + after-scan action = automated 1-D alignment, so you can easily implement an "Align" button.

## ...sscan module

---

- The sscan record
  - performs 1-D scan
  - before-scan link – optional completion callback
  - positioner: any writable, numeric, scalar PV (menus, enums are ok)
  - detector trigger: any writable, numeric, scalar PV
  - detector signal: any readable, numeric, scalar or 1D array PV
  - array detectors: exactly `<scanRecord>.NPTS` elements are acquired
  - array trigger: callback indicates array data are ready to read
  - after-scan link – optional completion callback
  - pause/resume
  - abort (`<scanRecord>.EXSC -> 0`) wait for callbacks, cleanup
  - kill (two aborts in a row) abandon callbacks
  - handshake with multiple display / data-acquisition clients
  - handshake with data-storage client

# Other data-acquisition-related software

---

- **Data-visualization tools for use with synApps**
  - Run-time look at scan data
  - Offline tools for data-file manipulation
  - Supports 1-3 dimensional data
  - Distributed independently of ioc software
  - See lecture “*Data Visualization.*”
  
- **CCD data-acquisition tools**
  - 1) CCD module (see lecture “*Detectors and Feedback*”)
  - 2) Portable CA Server based CCD support, and related software
    - <http://www.aps.anl.gov/aod/bcda/dataAcq/index.php>
  - Both of these solutions allow an EPICS CA client to drive data acquisition.
  - Both support `ca_put_callback()`, as required by the **sscan** record.



# std module

---

- **Epid record**
  - Extended PID record – see “*Detectors and feedback*” lecture
- **Scaler record**
  - Controls a set of counters with a common clock, gate, and trigger
- **String-sequence record**
  - Like the seq record in base, but works for strings and numbers
  - Can choose to wait for completion after each step in sequence
- **Soft-motor database (Jonathan Lang)**
  - *Run-time programmable* soft-motor/transform/hard-motor database
  - Quick solution for driving a motor through a nonlinear transform
- **Timestamp record [SLAC]**
  - needed by SNS’ vxStats; currently not available in a module
- **4-step database**
  - Up to four steps of (set condition; read data) with an end calculation
  - Originally developed for dichroism experiments

# utils directory

---

- **changePrefix**
  - Global search and replace of EPICS PV prefix within a copy of the xxx module
- **copyAdl**
  - Find all MEDM-display files buried in a file tree; copy to specified directory.

# vme module

---

- **VME record**
  - Provides run-time access to VME bus
  - Great for testing hardware
  - Run-time programmed control of an unsupported VME board
  
- **Device support for VME hardware**
  - Joerger scaler
  - APS bunch-clock generator
  - APS machine-status interface
  - Heidenhain encoder interpolator
  - Generic A32 VME interface
  - HP Laser interferometer
  - VMI4116 16-bit DAC
  - Acromag 9440 16-bit digital input

# xxx module

---

- **Prototype user directory**
  - Builds everything in synApps into a load module
  - Contains command files to load/configure everything in synApps
  - Contains sample top-level MEDM-display file
  - Contains sample script to set environment variables and start up the sample user interface
  - Contains table of recommended address/interrupt configuration for selected VME and IndustryPack hardware
  
- **Two ways to use this module**
  - 1) Make copies; run changePrefix; build; customize; run a beamline
    - *this is the recommended use*
    - *detailed instructions in support/documentation*
  - 2) Reference/grab bag

# For developers: features of synApps

---

- **extended-processing records**
  - records that are neither synchronous nor asynchronous, as these terms are described in the *EPICS Application Developer's Guide*
- **completion reporting**
  - All databases behave correctly when written to by `ca_put_callback()`.
- **recDynLink links**
  - Similar to standard EPICS links, but no “PP NMS” attributes
- **GUI standards**
  - Default colors for menus, PV values, links, etc.
- **coordinated motions**
  - Many of the databases in synApps (especially in ‘optics’) involve coordinated motion of several motors.
- **initialization of complex databases**
  - Some common EPICS initialization problems are handled in various synApps databases.

# Coordinated motions

---

- **Simple cases: database (transform records)**
  - Slits, mirrors, spherical-grating monochromator
- **More complicated cases: SNL code**
  - Multiple-crystal monochromators
- **Very complicated cases: custom record**
  - Optical table, scan
- **Criteria a useful coordination should meet:**
  - Report completion to `ca_put_callback()`
  - Share control of base positioners with CA clients
  - Recover state from the states of base positioners

# Completion reporting

---

- **Simple prescription for databases contained within a single ioc:**
  - Use only PP links and forward links in execution chain.
  
- **Database operations spanning more than one ioc:**
  - Use records with put\_callback links to span iocs:
    - *calcout* with asynchronous device support
    - *sscan, swait*
    - *sseq* or *sCalcout* (with *.WAIT\** = "Wait")
  
- **Cases in which a CA client performs part of the operation:**
  - 1) Database sets a **busy** record via PP or put\_callback link.
  - 2) CA client clears the **busy** record when operation is done.
  
- **Cases in which part of the operation is driven by a CP link:**
  - Not different from above; a CP link is a CA client

# Initialization of complex databases

---

- **Initial values: .VAL vs. .DOL**
  - Most records allow .VAL field to be set in the database.
  - Note that .DOL cannot be used for constant strings.
  
- **Save-restore and interaction with record/device initialization**
  - 1) save-restore pass 0
  - 2) record/device initialization → *device support can use pass-0 value*
  - 3) save-restore pass 1 → *pass-1 overrides record/device-init value*
  
- **.PINI (Process at INIt) uses and limitations**
  - This is the normal mechanism for database initialization.
  - What if you need a value from some other .PINI-initialized record, and that record hasn't processed yet?
  - Note .PHAS is not considered in .PINI processing.



# ...Initialization of complex databases

---

- **Contending with link alarms**
  - If you have an input link to a record with `.UDF=1`, you get a link alarm.
  - `.UDF=1` until a record processes. (In 3.14.1+, database can specify `.UDF`)
    - The transform record can abort execution on a link alarm (or not).
- **Initialization problems with CP links**
  - You have a CP link to a field that is a calculation result.
  - If the calc result is the same as the field's initial value, you'll have the right value, but you won't *know* that you have the right value, and you won't know for how long to wait to be sure.
    - The transform record *always* posts its initial calculation result.
- **Programmatically initializing link fields**
  - Link field must be written with a CA link (because lock-set recalculates).
  - `.PINI` processing occurs *before* CA is running (EPICS 3.13.5+).
  - Can't use `.PINI`; Drive init from a scan task; set init record to "Passive" when init is done.