




Getting Started with EPICS Applications / Special Topics

Introduction to synApps (v5.1)

Tim Mooney
1/11/2005

Argonne National Laboratory

A U.S. Department of Energy
Office of Science Laboratory
Operated by The University of Chicago



What is synApps?

- A collection of EPICS applications for synchrotron-beamline users
<http://www.aps.anl.gov/aod/bcda/synApps>
- EPICS modules and build/configuration tools:
 - Modules: *autosave, calc, camac, ccd, dac128V, dxp, ip, ip330, ipUnidig, love, mca, motor, optics, quadEM, sscan, std, vme, xxx*
 - Build/config: *config and utils directories*
- Related clients, libraries, and visualization tools:
 - IDL: *scanSee, mca display, ezcalDL, ezcaScan, ez_fit, HDF translator/browser, Ascii-format plotter, image processors, etc.*
 - CA-Server based CCD control
 - some python support

Pioneering Science and Technology

Getting Started with EPICS IOCs: Record Types and Examples

Office of Science U.S. Department of Energy



synApps modules

- Modules contain the following kinds of support:
 - Compiled code; libraries
 - *E.g., record and device support*
 - *State-Notation-Language programs*
 - EPICS databases and autosave-request files
 - *A database is a **program** written in a high-level language.*
 - *One or more copies of a database can be run, each with its own private variables (PV's).*
 - *The database designer recommends PV's to be autosaved by naming them in a .req file; you can override with a private copy of the file.*
 - MEDM-display files
 - *The default user interface*
 - Documentation

Pioneering Science and Technology

Getting Started with EPICS IOCs: Record Types and Examples

Office of Science U.S. Department of Energy


Other EPICS modules used by synApps

- **asyn 4.1**
- **ipac 2.8**
- **seq 2.0.8 (9?)**
- **genSub 1.6**
- **vxStats 1.7.2c**
- **allenBradley 2.1**

Pioneering Science and Technology

Getting Started with EPICS IOCs: Record Types and Examples

Office of Science U.S. Department of Energy



autosave module



- **Records latest values of selected EPICS PVs; restores those values when the ioc restarts.**
 - not an archiver; only the latest value is saved
 - not the same as saveData, which writes scan data
 - When a list of PV's is saved, the entire list is written, even if only one PV has changed.
- **Can save/restore any scalar or array-valued PV (synApps 5.1)**
 - Array-valued PV must be hosted by the ioc that does the restore operation. (Typically, all ioc's save/restore their own PV's.)
 - DBF_MENU, DBF_ENUM PV's are handled by number.
- **Save operation uses channel access for scalars.**
- **Restore operation uses static database access for scalars.**
- **Arrays are saved and restored with database access.**



...autosave module



- **Three restore options for save files:**
 - 1) before record/device initialization
 - Motor positions must be restored at this time.
 - Arrays cannot be restored at this time. *
 - PV's that are DBF_NOACCESS before record init (e.g., genSub variable-type fields) cannot be restored at this time. *
 - 2) after record/device initialization
 - to override record-initialization values
 - Link fields cannot be restored at this time. *
 - 3) both before and after record initialization
 - The 'auto_settings.sav' file is restored at both times.
 - It's not an error to attempt to restore a PV at the wrong time.
 - If you restore a motor position at this time, you override the value read from hardware, without writing to hardware.

* Not illegal, just doesn't work



...autosave module



- **PV lists can use include files (e.g., <database_name>.req), include path.**
 - Database developer can supply default include file with database.
 - User can override with custom include file.
- **Save triggers:**
 - on change of any PV in the list
 - periodically
 - on change of a trigger PV
 - manual
- **User can reload save sets.**
- **Autosave can recover from file-server reboot (synApps 4.6+).**
 - Currently, only on vxWorks
- **User can choose to save redundant files (synApps 4.6+).**
- **Autosave reports status via EPICS PV's (synApps 5.1+).**



...autosave module



- **Sample request file**

```
xxx:my_PV.VAL
xxx:my_array_PV.VAL
file motor_settings.req P=${P},M=m1
...
<END><lf>
```

PV names
Keyword
Macro replacements
Name of include file

- **Sample save file**

```
# save/restore V4.4 Automatically generated..
xxx:my_PV.VAL 1.0
xxx:my_array_PV.VAL @array@ { "0" "0.1" ... "10.2" }
xxx:m1.DIR 0
xxx:m1.DHLM 100
xxx:m1.DLLM -100
...
```



calc module

- Support for evaluation of string or numeric expressions entered at run time (or at database-configure time)
- Records
 - **sCalcout** – like calcout, but also supports string expressions; user can specify wait-for-completion.
 - **swait** – like calcout, but uses recDynLink (no “PP MS” link attributes)
 - **transform** – like 16 calcout records that share a PV data pool
- Other code
 - string-calc engine
 - sCalcout soft device support (with wait-for-completion option)
 - interpolation routines for genSub record
 - (yet another) averaging routine for sub record



...calc module

- Databases, medm displays for run-time programming
 - userCalc
 - userStringCalc
 - userTransform
 - userAve
 - arrayTest
 - interpolation
- Examples of ALL calc expressions (normal and stringCalc) can be found in synApps MEDM help displays



camac module

- Communication with CAMAC crate/modules
- Records
 - camac – generic BCNAF/data for run-time camac control
- Devices supported
 - VME bus adapter
 - CAMAC crate controller
 - E500 motor controller
 - RTC-018 real-time clock
 - QS-450 quad scaler
 - DXP spectroscopy system (now in *dxp* module)



ccd module

- Support for area detectors (CCD's and image plates)
- Supported devices
 - MAR 165 CCD
 - MAR 345 image-plate reader
 - Roper (all WinView-supported CCD's, including former Princeton and most former Photometrics devices)
 - Bruker SMART CCD
- Can control, at minimum
 - exposure time
 - file name
 - data-acquisition start
 - wait for acquisition to complete
 - much more for most devices
- See lecture “*Detectors and Feedback.*”



config directory



- **Configures and builds all modules in or used by synApps**
- **MASTER_RELEASE**
 - specifies version number and file path to EPICS base, and to every module in or used by synApps
- **makeReleaseConsistent.pl**
 - Edits <module>/configure/RELEASE for every module in or used by synApps, to agree with MASTER_RELEASE
 - “gnumake release” causes this to run.
- **Makefile**
 - “gnumake <whatever>”, in config directory, does <whatever> for all modules.



Pioneering
Science and
Technology

Getting Started with EPICS IOCs: Record Types and Examples



Office of Science
U.S. Department
of Energy

dac128V module



- **device support, database, and MEDM displays for dac128V IndustryPack module**
 - 8-channel, 12-bit DAC
 - Support exists to run a DAC channel manually, or according to an algorithm written at run time, or as a *scan* positioner, or as part of a PID feedback loop.
- **See lecture “*Detectors and Feedback.*”**



Pioneering
Science and
Technology

Getting Started with EPICS IOCs: Record Types and Examples



Office of Science
U.S. Department
of Energy

documentation directory



- **TOP-level synApps documentation**
 - What synApps is
 - How to build it
 - How to make a user application from the ‘xxx’ sample module
 - How to fit the user application to a particular set of hardware
- **This presentation**



Pioneering
Science and
Technology

Getting Started with EPICS IOCs: Record Types and Examples



Office of Science
U.S. Department
of Energy

dxp module



- **record, device support, databases, and MEDM displays for XIA DXP and Saturn spectroscopy systems**
- **dxp record for setting DXP parameters**
- **device support for the mca record**
- **See lecture “*Detectors and Feedback.*”**



Pioneering
Science and
Technology

Getting Started with EPICS IOCs: Record Types and Examples



Office of Science
U.S. Department
of Energy

ip module

- **device support, SNL code, databases, and MEDM displays for many message-based devices**
 - originally, for devices supported via IndustryPack hardware
 - Note some of this support will inevitably be out of date -- pending access to hardware for testing.
- **deviceCmdReply (was serial_OI_block, GPIB_OI_block)**
 - Used to write support at run time for one command/reply message
 - sCalcout to format output string
 - asyn record to write/read device
 - sCalcout record to parse reply
- **devXxStrParm device support**
 - probably will be replaced by streams/asyn

ip330 module

- **device support, databases, and MEDM displays for the IP330 ADC IndustryPack module**
- **16/32 channel, 16-bit ADC**
 - ip330Scan for periodic, averaged reads of ADC channels
 - ip330Sweep, with the MCA record, for using ip330 as a waveform-digitizer
 - ip330PID for using the ip330 in a fast-feedback loop
- **See lecture “Detectors and Feedback.”**

ipUnidig module

- **device support, databases, and MEDM displays for the IPUnidig digital I/O IndustryPack module**
- **IP-UD-I 24-channel input/output/interrupt module**
- **DIO316I 48-bit digital I/O module**
- **See lecture “Detectors and Feedback.”**

love module

- **Support for Love controllers**
 - orphaned, currently under re-development for EPICS 3.14

mca module

- **Support for multichannel analyzers, multichannel scalers, and other array-valued detectors**
- **mca record**
- **device support**
 - Canberra 556 AIM module (MCA and ICB controller)
 - DSA-2000 Ethernet MCA
 - various Canberra-ICB modules for spectroscopy
 - SIS 3801 (Struck STR7201) MCS
 - (DXP support in dxp module)
 - (IP330 support in ip330 module)
 - (quadEM support in quadEM module)
- **See lecture “Detectors and Feedback.”**

motor module

- **Motor record and device support**
 - stepper and servo motors
 - soft-motor support
 - *Put motor “face” on, e.g., a DAC channel*
 - *Drive a hard motor through a nonlinear transform*
 - user/dial/raw coordinates
 - backlash-takeout algorithm
 - pre/post move commands
 - many more features
- **See lecture “Motors.”**

optics module

- **Slits and mirrors**
 - *Four virtual positioners; two real motors*
 - *Automatic sync to motor positions*
 - *Completion reporting*
- **Monochromators**
 - *Nondispersive double-crystal*
 - Geometries: (Y1, Z2), (Y2, Z2)
 - Crystal species: Si, Ge, Diamond, Si (77K)
 - Miller indices, allowed reflections
 - Operational modes:
 - *Use/Set*
 - *Manual/Auto*
 - Managing the vertical beam offset
 - Automatic sync to motor positions

...optics module

- **...Monochromators**
 - *Spherical grating*
 - Geometrical variables:
 - 1) *Grating line density; radius*
 - 2) *Tangent-arm length*
 - 3) *Diffraction order*
 - 4) *Input/output slit distances*
 - Operational modes:
 - *Use/Set*
 - *Manual/Auto*
 - Grating-stripe list
 - Manual sync to motor positions

...optics module



- **...Monochromators**
 - *Dispersive double-crystal*
 - Geometries: nested, symmetric
 - Crystal species: Si, Ge, Diamond, Si (77K)
 - Miller indices, allowed reflections
 - Operational modes:
 - *Use/Set*
 - *Manual/Auto*
 - *Theta 1 / Theta 1&2 / Rock Theta2*
 - Accommodate incident-beam angle shift (“world offset”)
 - Automatic sync to motor positions

...optics module




- **Optical table**
 - **Table** record supports a six-degree-of-freedom optical table.
 - User/client can write either to $(x, y, z, \theta_x, \theta_y, \theta_z)$, or to underlying motor records.
 - Table rotates about user-specified point.
 - Table database includes a list of rotation points, selected by menu.
 - Can recover table position from motor positions
 - Partial support for fewer than six degrees of freedom

...optics module



- **Optical table**
 - Geometries
 - SRI
 - GeoCARS
 - Newport
 - PNC
- Calibration/sync
 - Use/Set – changes to $[X, Y, ..]$. move table / change calibration
 - Zero – redefine current $[X, Y, ...]$ as zero
 - Sync – update $[X, Y, ...]$ from motors, honoring calibration
 - Init – clear calibration and sync to motors
- Table record sets motor speeds so that motors start/stop together.



quadEM module



- **Support for APS Detector Group's (Steve Ross) four-input electrometer.**
- **See lecture “Detectors and Feedback.”**

sscan module



- **Support for user-programmable data-acquisition**
 - **sscan** and **busy** records
 - saveData
 - recDynLink
- **A one-dimensional scan:**
 - Do NPTS times:
 - Set conditions e.g., move motors; wait for completion
 - Trigger detectors e.g., start scaler; wait for completion
 - Acquire data read detector signals; store in arrays
 - Write data to NFS file
- **Multidimensional scan:**
 - Same as a 1-D, but detector trigger executes inner-loop scan.
 - saveData monitors a set of **sscan** records, determines scan dimension when scan starts, and writes data as it is acquired.



Getting Started with EPICS IOCs: Record Types and Examples



...sscan module



- **scan features:**
 - Three 1-D scan types: constant-step-size, table-driven, fly
 - Unlimited number of data points, scan dimensions
 - 0-4 positioners, 0-4 detector triggers, 0-70 detector signals
 - Acquisition from scalar and 1-D-array-valued PV's
 - Detector/client wait, data-storage wait
 - Pause/resume, abort
 - Double buffered: can write 1-D acquired data during next 1-D scan
 - saveData writes self-describing XDR-format (".mda") files to NFS-mounted disk (vxWorks only, at present).
 - A positioner can have private scan parameters (scanparm record).
 - After-scan actions include move to peak, valley, and edge.
 - scanparm record + after-scan action = automated 1-D alignment, so you can easily implement an "Align" button.



Getting Started with EPICS IOCs: Record Types and Examples



...sscan module



- **The sscan record**
 - performs 1-D scan
 - before-scan link – optional completion callback
 - positioner: any writable, numeric, scalar PV (menus, enums are ok)
 - detector trigger: any writable, numeric, scalar PV
 - detector signal: any readable, numeric, scalar or 1D array PV
 - array detectors: exactly `<scanRecord>.NPTS` elements are acquired
 - array trigger: callback indicates array data are ready to read
 - after-scan link – optional completion callback
 - pause/resume
 - abort (`<scanRecord>.EXSC -> 0`) wait for callbacks, cleanup
 - kill (two aborts in a row) abandon callbacks
 - handshake with multiple display / data-acquisition clients
 - handshake with data-storage client



Getting Started with EPICS IOCs: Record Types and Examples



Other data-acquisition-related software



- **Data-visualization tools for use with synApps**
 - Run-time look at scan data
 - Offline tools for data-file manipulation
 - Supports 1-3 dimensional data
 - Distributed independently of ioc software
 - See lecture "**Data Visualization.**"
- **CCD data-acquisition tools**
 - 1) CCD module (see lecture "**Detectors and Feedback**")
 - 2) Portable CA Server based CCD support, and related software
 - <http://www.aps.anl.gov/aod/bcda/dataAcq/index.php>
 - Both of these solutions allow an EPICS CA client to drive data acquisition.
 - Both support `ca_put_callback()`, as required by the **sscan** record.



Getting Started with EPICS IOCs: Record Types and Examples



std module

- **Epid record**
 - Extended PID record – see “*Detectors and feedback*” lecture
- **Scaler record**
 - Controls a set of counters with a common clock, gate, and trigger
- **String-sequence record**
 - Like the seq record in base, but works for strings and numbers
 - Can choose to wait for completion after each step in sequence
- **Soft-motor database (Jonathan Lang)**
 - *Run-time programmable* soft-motor/transform/hard-motor database
 - Quick solution for driving a motor through a nonlinear transform
- **Timestamp record [SLAC]**
 - needed by SNS’ vxStats; currently not available in a module
- **4-step database**
 - Up to four steps of (set condition; read data) with an end calculation
 - Originally developed for dichroism experiments



utils directory

- **changePrefix**
 - Global search and replace of EPICS PV prefix within a copy of the xxx module
- **copyAdl**
 - Find all MEDM-display files buried in a file tree; copy to specified directory.



vme module

- **VME record**
 - Provides run-time access to VME bus
 - Great for testing hardware
 - Run-time programmed control of an unsupported VME board
- **Device support for VME hardware**
 - Joerger scaler
 - APS bunch-clock generator
 - APS machine-status interface
 - Heidenhain encoder interpolator
 - Generic A32 VME interface
 - HP Laser interferometer
 - VMI4116 16-bit DAC
 - Acromag 9440 16-bit digital input



xxx module

- **Prototype user directory**
 - Builds everything in synApps into a load module
 - Contains command files to load/configure everything in synApps
 - Contains sample top-level MEDM-display file
 - Contains sample script to set environment variables and start up the sample user interface
 - Contains table of recommended address/interrupt configuration for selected VME and IndustryPack hardware
- **Two ways to use this module**
 - 1) Make copies; run changePrefix; build; customize; run a beamline
 - *this is the recommended use*
 - *detailed instructions in support/documentation*
 - 2) Reference/grab bag



For developers: features of synApps



- **extended-processing records**
 - records that are neither synchronous nor asynchronous, as these terms are described in the *EPICS Application Developer's Guide*
- **completion reporting**
 - All databases behave correctly when written to by `ca_put_callback()`.
- **recDynLink links**
 - Similar to standard EPICS links, but no "PP NMS" attributes
- **GUI standards**
 - Default colors for menus, PV values, links, etc.
- **coordinated motions**
 - Many of the databases in synApps (especially in 'optics') involve coordinated motion of several motors.
- **initialization of complex databases**
 - Some common EPICS initialization problems are handled in various synApps databases.



Coordinated motions



- **Simple cases: database (transform records)**
 - Slits, mirrors, spherical-grating monochromator
- **More complicated cases: SNL code**
 - Multiple-crystal monochromators
- **Very complicated cases: custom record**
 - Optical table, scan
- **Criteria a useful coordination should meet:**
 - Report completion to `ca_put_callback()`
 - Share control of base positioners with CA clients
 - Recover state from the states of base positioners



Completion reporting



- **Simple prescription for databases contained within a single ioc:**
 - Use only PP links and forward links in execution chain.
- **Database operations spanning more than one ioc:**
 - Use records with `put_callback` links to span iocs:
 - *calcout* with asynchronous device support
 - *sscan*, *swait*
 - *sseq* or *sCalcout* (with `.WAIT*` = "Wait")
- **Cases in which a CA client performs part of the operation:**
 - 1) Database sets a **busy** record via PP or `put_callback` link.
 - 2) CA client clears the **busy** record when operation is done.
- **Cases in which part of the operation is driven by a CP link:**
 - Not different from above; a CP link is a CA client



Initialization of complex databases



- **Initial values: .VAL vs. .DOL**
 - Most records allow `.VAL` field to be set in the database.
 - Note that `.DOL` cannot be used for constant strings.
- **Save-restore and interaction with record/device initialization**
 - 1) save-restore pass 0
 - 2) record/device initialization → *device support can use pass-0 value*
 - 3) save-restore pass 1 → *pass-1 overrides record/device-init value*
- **.PINI (Process at INIt) uses and limitations**
 - This is the normal mechanism for database initialization.
 - What if you need a value from some other `.PINI`-initialized record, and that record hasn't processed yet?
 - Note `.PHAS` is not considered in `.PINI` processing.





...Initialization of complex databases

- **Contending with link alarms**
 - If you have an input link to a record with .UDF=1, you get a link alarm.
 - .UDF=1 until a record processes. (In 3.14.1+, database can specify .UDF)
 - The transform record can abort execution on a link alarm (or not).
- **Initialization problems with CP links**
 - You have a CP link to a field that is a calculation result.
 - If the calc result is the same as the field's initial value, you'll have the right value, but you won't *know* that you have the right value, and you won't know for how long to wait to be sure.
 - The transform record *always* posts its initial calculation result.
- **Programmatically initializing link fields**
 - Link field must be written with a CA link (because lock-set recalc).
 - .PINI processing occurs *before* CA is running (EPICS 3.13.5+).
 - Can't use .PINI; Drive init from a scan task; set init record to "Passive" when init is done.