


Data Visualization

EPICS Extensions IDL Tools

Ben-chin K. Cha
Beamline Control & Data Acquisition

Argonne National Laboratory



A U.S. Department of Energy
Office of Science Laboratory
Operated by The University of Chicago



Introduction

- Written in IDL/ezcaIDL to view real-time / post scan data
- Packaged data browsers with complete mouse driven user interface
- IDL interactive graphics and analysis tools
- High quality interactive graphic tools (IDL 6.X)
- Packaged tools support license free IDLVM 6.0
- Platforms: Unix , Win2000/XP, Linux
- Object methods for easy post scan data array extraction and manipulation



2


Outline

- Overview
- System and setup requirement
- ScanSee program user interface
- Sharable visualization subprograms
- Scan object methods and examples
- EzcaIDL setup and examples

3



Scan Visualization Architecture



```

graph TD
    ScanSee[ScanSee: Packaged IDL Scan Visualization Analysis Tools] <--> ezcaIDL[ezcaIDL]
    ezcaIDL <--> EzcaScan[EzcaScan/ezca]
    EzcaScan <--> IOC[IOC Controller sscan with saveData]
    IOC <--> EPICS[EPICS CA]
    EPICS <--> MEDM[MEDM Scan Controls]
    IOC --> Save[Save MDA Data]
    Save --> ScanSee
  
```

4

Overview

- **Scan Visualization Tools**
 - *scanSee* – 85 detectors, 1D/2D/3D scan (read MDA)
 - *catcher* – 15 detectors, 1D/2D scan (save/read 1D/2D files)
- **HDF Visualization Tools**
 - *hdfb* – HDF4 1D/2D/3D browser
 - *h5b* – HDF5 1D/2D/3D browser
- **Common/sharable Visual/Analysis Tools** <http://www.aps.anl.gov/~cha/dl.html/list.html>
 - *plot1d*, *ez_fit*, *overlay_1d*
 - *plot2d_image*, *plot2d*, *image2d*, *view3D_2D*
 - *calibration_factor*, *scan2d_roi*, *multiroi_pick*
 - *panimage*, *scan2d_overlay*, *wd_readascii*

Input Forms

- **XDR catcher 1D/2D file**
- **MDA 1D/2D/3D scan data**
- **Nexus/HDF release 4 1D/2D/3D data**
- **Nexus/HDF release 5 1D/2D/3D data**
- **Arrays data**
- **Fix format ASCII data**

Output Forms

- **1D/2D graphic output**
- **1D/2D ASCII report data**
- **Graphic output: TIFF / PNG / PICT / PS**
- **1D/2D ROI statistics reports**
- **Flexible 3D to 2D, 2D to 1D sliced output**

System Requirement

- **EPICS extensions built: ezca, EzcaScan, ezcalDL, idl**
 - *Scripts and executables installed under*
`/usr/local/epics/extensions/bin/solaris-sparc for 3.14.X`
 - *IDL programs installed under*
`/usr/local/epics/extensions/idllib for 3.14.X`
- **IDL 6.0 and IDLVM 6.0 installed**
- **IOC sscan record properly configured by medm**

Unix Setup Requirement

- Set EPICS_EXTENSIONS

```
setenv EPICS_EXTENSION /usr/local/epics/extensions
```
- Access EPICS 3.14.X

```
setenv EPICS_HOST_ARCH solaris-sparc
set path=( /usr/local/epics/extensions/bin/solaris-sparc $path)
```

Packaged Unix Scripts for IDL Visualization Tools

- **scanSee** – MDA viewer with real-time scanning
- **sscanscan** – post scan viewer for MDA files
- **catcher** – real-time scanning with data saving
- **viewer** – post scan viewer for data catcher
- **hdfb** – data browser for HDF4/Nexus files
- **h5b** – data browser for HDF5/Nexus files
- **img** – Tiff/Jpeg/Png/Ascii/Xdr image processor
- **idlvms** – run IDLVM 6.0 with any IDL 6.0 saved programs

IDLVM 6.0

- **No license required for IDL 6.0 saved program**
- **Document:** <http://www.aps.anl.gov/~chafid/vm60.html>
- **IDL 6.0 *.sav files:**
 - *catcher.sav* – for catcher
 - *viewer.sav* – for viewer without CA
 - *sscanscan.sav* – for scanSee R3.4 and later
 - *SB2.sav* – for MDA reader without CA
 - *hdfb.sav* – for HDF version 4
 - *h5b.sav* – for HDF version 5
 - *img.sav* – image processor
- **Unix script methods:**

```
idlvms <pname>
```

where <pname> can be any of the above name without *.sav
- **Tools & data calibration not available in IDLVM 6.0**


Unix Access Methods

- **IDLVM version** (recommend unless calibration is desired)

```
idlvms sscanscan          (scanSee.R3.4 and later)
idlvms catcher
```
- **Runtime version**

```
scanSee
catcher
```
- **Developer Version** (ITools available)

```
scanSee -D
catcher -D
```



Files Used by scanSee

- MDA scan files automatically saved by IOC
- Restart Configuration file: DC.config
- Saved output files:
 - ASCII/*.txt - 1D/2D data files
 - TIFF/*.tiff - TIFF image files
 - PNG/*.png - PNG image files
 - PICT/*.pict - MAC image files
 - ROI/*.roi, *.rpt - various ROI files
 - idl.ps, plot2d.ps ... - saved various PS plot, report files

Features Access Restriction


- ITOOLS
 - Developer version only
 - License free ITOOLS
 - 'readascii.sav' (IDLVM 6.1)
 - Idlvm readascii
- 1D/2D Calibration features
 - Developer and Runtime versions

Invocation Problem

- Problem of startup scanSee
 - Due to MDA file not found in configuration file
 - Wrong type or bad file picked
- Resolve Invocation Problem
 - Method 1
 - Use File->Open to pick a new file
 - Method 2
 - Remove the 'DC.config' before run scanSee


General Widget Interfaces

- **Click LMB** - selection of menu/button/cursor/slider/list
- **Text field** - text entry must be ended with CR
- **Button...** - a pop up window/subprogram associated with it
- **Multiple list selection**
 - *Cntl+LMB* – add/delete item to selection list
 - *Shift+LMB* – select all items between current and last click
- **ROI Mouse buttons**
 - *LMB* – pick/drag ROI
 - *MMB* – resize ROI
 - *RMB* – accept ROI



ScanSee Features

- **Use file DC.config** for easy restart
- **Flexible MDA file selection dialog**
- **User settable monitor of 1D/2D scan record**
- **Display near real-time 1D and 2D scanning data**
- **Flexible 1D display in detectors, line texture, plot ranges and labels, zoom, FWHM, easy fitting, etc**
- **Image data analysis, ROI statistics, calibrations**
- **Easy access 1D/2D/3D post scan data and pass to appropriate subprograms**
- **Easy ASCII/TIFF/PNG/PIC/ROI/PS output**



ViewData->
IMAGE2D..

Set New 2D Scan Ranges

Set New 2D Scan ranges (Yes/No)?

Xmin:	2idd:scan1.P1SP=-639.20000
Xmax:	2idd:scan1.P1EP=-625.30000
Ymin:	2idd:scan2.P1SP=-850.20000
Ymax:	2idd:scan2.P1EP=-834.30000

Yes No

21

Pioneering Science and Technology U.S. Department of Energy

PanImage Selection Dialog

Select Sublist: D01-D10 Grid Exp Factor: 1

Select Detectors: D01 D02 D03 D04 D05

DI/Row: 10

Accept All Color... Close

Save PanImage Type TIFF Reverse ASCII...

Filename: pan.tiff

PanImage_sel Features

- Panimage Select sub-list: D01-D10,D11-D20,...
- Multiple detectors selection list, default All
- Grid factor, default 1 (60x60)
- Show/Hide detector names
- Detectors per row, default 10
- Accept, All, Color table dialog
- Easy TIFF, PNG, PICT, ASCII panimage generation

22

Pioneering Science and Technology U.S. Department of Energy

Panimage_sel Output

Scan # 2

Color... Printer... Print PS Scale 1 Done

23

Pioneering Science and Technology U.S. Department of Energy

2D Calibration_factor Features

Calibration - Scan # 2

File Help Done

Calibration Factor

F1	1.00000
F2	1.00000
F3	1.00000
F4	1.00000
F5	1.00000
F6	1.00000
F7	1.00000
F8	1.00000
F9	1.00000
FA	1.00000
FB	1.00000
FC	1.00000
FD	1.00000
FE	1.00000
FF	1.00000
FDL	1.00000


Default: Accept ReCalc... ASCII... PICT... PICT...

of Over: 4 F1409-D7

24

Pioneering Science and Technology U.S. Department of Energy

2D Calibra Result




Plot2d Features


- Support TV, Surface, Contour, Shade Surface
- Access of ROI, Pick1D, ASCII data sub-programs
- TV image mouse query and sub-region query
- Dialog for various plot area configurations
- Dialog for various TV manipulation features: sharp, smooth, hist_equal,

Plot2d_TV Options

- TV-log, TV, TVSCL
- Image query mode
- Hist_Equal
- Max/min filter value
- Smooth/Unsharp Mask
- Roberts/Sobel Sharpping
- Threshold




Scan2d_roi



ROI Features

- Rect/Filter/Polygon/Multi mode of ROI
- Various controls of different ROI mode
- Mouse query / modify ROI
- Display statistic summary results
- User Settable ROI define and report filename
- Dialogs save, display, manage ROI reports


MultiROI_pick subprogram



Multi_ROI Features


- Flexible Add/delete/modify ROIs
- Free hand drawing of polygon ROIs
- Display one/all statistic of ROIs
- Save/load multi ROIs to/from a file
- Manipulate/generate multi ROIs statistic report
- Save ROIs as a TIFF/PNG/PICT file
- Image query

ViewData->SSCAN (1D/2D/3D/4D)...




View3d_2d Slicer

Features



- Plot2d/Ascii2d/ROI2d/Pick1d subprogram selection
- Flexible axial rank selection for 2D slicer
- Flexible 2D image animator or preview slider
- 2D image slice # selection
- View3d_2d Sum ROI subprogram



View3d_Sum_2DROI Features

- Slider controls for I,J, axial range K1,K2
- Spectrum plot at I,J with K1,K2 marked as blue lines
- Option of PLOT1D, PLOT1D_ALL_I, PLOT1D_ALL_J spectrum
- Image of sum of 2D values bound by K1,K2
- Plot of cross image line values at I,J
- Option of Plot2d/Ascii2d/Roi2d/Pick1d/Data/PS plot

CALIBRA_PICK1D_2idd_0002.rmda : 3D Seq # :D_16X_Slice.slc0_19


	Col27	Col28	Col29
Row0	1111,00	1190,00	1344,00
Row1	1313,00	1412,00	1617,00
Row2	2233,00	2532,00	3151,00
Row3	3205,00	3309,00	4062,00
Row4	5700,00	5882,00	6988,00
Row5	7558,00	7450,00	7677,00
Row6	9674,00	8316,00	8738,00
Row7	9655,00	9533,00	9777,00
Row8	10631,0	10897,0	10619,0
Row9	11304,0	11473,0	11039,0

Calibra_pick1d Dialog

Pick1D Features

- Table display of 2D image value with format control
- Column/Row selection pops up plot1d program
- Plot all columns/rows with plot1d subprogram
- Access ez_fit subprogram with image array

Plot All I-spectrum



Plot1d Features

- Plot region resize with window size
- Export ASCII, TIFF, PNG, PICT, PS file generation
- Access of IPLOT of IDL 6.0 Itools (developer version)
- Multiple lines plot program support scattering data
- Options... Dialog for various line plot configurations, legend, statistic, fitting, FWHM, etc

EZ_FIT (R1.6)

File GetData Curve Fit Multi Fit Help Close

1D XRF File Fitting,bin

Open Vector CON
Open Vector CUR
Print Vector ERR
Close Vector Z GAU
2D Image Vector Z LA
POL
LORE
SV

POLY_FIT

Fitting, rpt... Rvs Bkg Options... Export Print IPLOT... Close

Least-Square POLY_FIT

Y = A0 + A1 * X + A2 * X^2 + A3 * X^3 + A4 * X^4

A0=0.00000000
A1=1.1941514
A2=-4.5051098
A3=0.7055048
A4=0.006784001
SQRM=782.5014

GOODNESS OF FIT = 745.39918

User Name: cha

EZ_FIT Dialog

Plot1d

Rvs Bkg Options... Export Print IPLOT... Close

Rank=0, I=27, J=27

Options...->FITTING...
Muti Fit->LORENTZIAN

LORENTZIAN MULTI-ROI Fitting


Define Multi Lorentzian ROIs / Curve Fitting

ROI #	Left	Right	Peak Intensity	Centroid # X	FWHM
5	0,000	1999	30,999225	367,92469	26,4492
0		1999		0	


Subplot	Y	ROI	XL	XR	Intensity	Centroid	FWHM/2
0	0,00000	1999,00	139,999	267,223	11,2246		
1	0,00000	1999,00	164,445	327,963	11,2246		
2	0,00000	1999,00	35,5317	363,353	11,2246		
3	0,00000	1999,00	55,6539	767,770	11,2246		
4	0,00000	1999,00	39,2629	945,640	11,2246		
5	0,00000	1999,00	39,3963	967,634	11,2246		

- Click First Peak
- Click Add Btn
- Click Next Peak
- Click Show ROIs
- Click Calc Fits...

Calc Fits...



idlvm readascii



Use Sharable Programs

- **IDL Developer License Required**
- **Common Visualization Tools Sharable by Any Other IDL Program**
 - *Plot1d, plot2d, image2d, view3d_2d, panimage,...*
- **IDL Startup Setup**
setenv IDL_STARTUP /usr/local/epics/extensions/idllib/viewer_startup.pro
- **Example**

```
idl
IDL> loadct, 39
IDL> da = dist(100)
IDL> plot2d, da
IDL> plot1d, da
```

Scan Object Files

- **Scan Object files with methods:**
 - *sscan__define.pro* (access MDA data)
 - *scan2d__define.pro* (access catcher 2D image file)
 - *scan1d__define.pro* (access catcher 1D file)
- **Object Access Setup :**
setenv IDL_STARTUP /usr/local/epics/extensions/idllib/viewer_startup.pro

Sscan Object Examples

- **Referring**
http://www.aps.anl.gov/~cha/idl_html/sscan_ObjectRef.html
- **Examples:**
idl
IDL> file='/home/beams/CHA/Yorick/data/2ixm_0020.mda'
IDL> V = obj_new('sscan',file=file)
IDL> V->image2d ; use object method
IDL> V->view3d_2d ; use object method
IDL> V->read,da1d=da1d,da2d=da2d,da3d=da3d,...
IDL> view3d_2d,da3d ; call sharable program
IDL> image2d,da2d ; call sharable program

Catcher: Scan1d Object Examples

- **Referring:**
http://www.aps.anl.gov/~cha/idl_html/scan1d_ObjectRef.html
- **Examples:**
idl
IDL> file='/home/beams/CHA/Yorick/data/junk2')
IDL> .run scan1d__define
IDL> V = obj_new('scan1d',file=file)
IDL> V->read,10,pa=pa,da=da,/plot; read scan #10 in pa and
; da array also plot da
IDL> V->plot,10,iy='4,5' ; read scan # 10 and plot
; detector 4 & 5 only
IDL> plot1d,da ; call sharable program

Catcher: Scan2d Object Examples

- **Referring:**
http://www.aps.anl.gov/~cha/idl_html/scan2d_ObjectRef.html
- **Examples:**
idl
IDL> file='/home/beams/CHA/Yorick/data/junk2.image')
IDL> V = obj_new('scan2d',file=file)
IDL> V->print ; summary of image file
IDL> v->panimage,2 ; access the 2nd scan
IDL> v->images,2,im,def ; extract all images in 2nd scan
IDL> image2d,im,id_def=def ; use image2d to view 2nd scan
; images

HDF4: NX Object Examples

- **Referring:**
http://www.aps.anl.gov/~cha/idl_html/NX_ObjectRef.html
- **Examples:**
idl
IDL> file='/home/beams/CHA/Yorick/data/2xfm_0020.nexus')
IDL> .run NX__define
IDL> v = obj_new('NX',file=file)
IDL> v->print ; dialog show file info, find num SDS
IDL> for i=0,309 do begin v->sds,i,da,name=na & help,da,out=t & print,na,i,t & end
; get all SDS data name and array info
IDL> retail
IDL> v->sds,131,da,name=na ; read the 131th set of SDS
IDL> view3d_2d,da ; pass 3D array to view3d_2d program
IDL> v->sds,96,da,name=na ; read the 96th set of SDS
IDL> plot2d,da ; pass 2D array to plot2d program
IDL> hdfb ; call hdfb program

Access ezcaIDL functions

- **Reference Documents:**

<http://www.aps.anl.gov/~cha/catcher/ezcaIDLRef.html>

- **Unix Setup Requirement:**

```
source /usr/local/epics/extensions/bin/solsris-sparc/ezcaidl_setup  
( LD_LIBRARY_PATH, EZCA_IDL_SHARE, IDL_STARTUP)
```

- **Examples:**

Idl -32

```
IDL> r = caGetArray('cha:scan1.D1DA',data,max=11)           ; get 11 data points  
IDL> plot1D,data                                           ; plot data array  
IDL> x = ['chademoai1', 'chademoai2']                     ; set x to two PVs  
IDL> r = caGetArray(x, y) & print,y                       ; get & print values  
IDL> newy = y+1. & r = caPutArray(x, newy)                 ; add & put values
```

45

Where to Get More Information

- <http://www.aps.anl.gov/aod/bcda/dataVis/index.php>

46